## INTENDED USE

### FOR IN VITRO DIAGNOSTIC USE

CLINIQA® Liquid QC™ Cardiac Marker Control (CK-MB/Myoglobin/Troponin I) is intended for use as an assayed qualitative quality control material for CK-MB, Myoglobin and Troponin I as well as a quantitative quality control material for Troponin I. CLINIQA Liquid QC Cardiac Marker Control (CK-MB/Myoglobin/Troponin I) is not intended for use as a standard.

### SUMMARY

In clinical laboratories the use of quality control materials to objectively monitor the accuracy and precision of procedures is well established. CLINIQA Liquid QC Cardiac Marker Control (CK-MB/Myoglobin/Troponin I) is provided at two levels to assist in the monitoring of analytical systems within the clinical range.

### PRODUCT DESCRIPTION

CLINIQA Liquid QC Cardiac Marker Control (CK-MB/Myoglobin/Troponin I) is prepared from human plasma, human serum and human proteins. Preservatives, stabilizers and sodium azide have been added to maintain product integrity. CLINIQA Liquid QC Cardiac Marker Control (CK-MB/Myoglobin/Troponin I) is a ready-to-use liquid control requiring no reconstitution or frozen storage.

### STORAGE AND STABILITY

CLINIQA Liquid QC Cardiac Marker Control (CK-MB/Myoglobin/Troponin I) is stable until the expiration date on the vial label when stored unopened at 2-8°C. Once opened, CLINIQA Liquid QC Cardiac Marker Control (CK-MB/Myoglobin/Troponin I) is stable for 30 days when stored tightly capped at 2-8°C.

### PROCEDURE

CLINIQA Liquid QC Cardiac Marker Control (CK-MB/Myoglobin/Troponin I) should be treated the same as patient specimens and run in accordance with the instructions accompanying the test system being used. Gently mix the contents of each vial before sampling to ensure homogeneity. Replace cap immediately and store at 2-8°C.

### LIMITATIONS OF PROCEDURE

CLINIQA Liquid QC Cardiac Marker Control (CK-MB/Myoglobin/Troponin I) should not be used past the expiration date on the vial label. If there is evidence of microbial contamination or excessive turbidity in the product, discard the vial.

CLINIQA Liquid QC Cardiac Marker Control (CK-MB/Myoglobin/Troponin I) is a stabilized

liquid product. To obtain consistent assay values, CLINIQA Liquid QC Cardiac Marker Control (CK-MB/Myoglobin/Troponin I) requires storage and handling as detailed in STORAGE and STABILITY.

Accurate and reproducible results are dependent upon properly functioning instruments and reagents. The published assay values were obtained using reagents and procedures available at the time of assay. In the event reagents or procedures are altered by the manufacturer, different assay values may be obtained.

### ASSIGNMENT OF VALUES

The assigned mean values were derived from analyses of vials representative of the entire lot. Analyte values were obtained from laboratories of the instrument manufacturer, manufacturers of instrument specific reagents, or from reference laboratories.

The Expected Range of the Mean is provided to assist the laboratory until it has established its own mean and standard deviation. It is considered good laboratory practice for each laboratory to establish its own mean and standard deviation for its test methods. The indicated Mean and Expected Range of the Mean should serve as a guide in assessing the performance of each test method.

The values are usually method dependent. The variations which can occur over time and between laboratories may be attributed to differences in laboratory technique, instrumentation, reagent lot, method modifications, and other systematic errors including random errors.

### CAUTION

#### Human source material. Treat as potentially infectious.

Each serum/plasma donor unit used in the manufacture of this product has been tested by FDA accepted methods and found non-reactive for the presence of HBsAg and antibody to HIV-1/2, HCV and HIV-1 Ag. While these methods are highly accurate, they do not guarantee that all infected units will be detected. Because no known test method can offer complete assurance the hepatitis B virus, hepatitis C virus, human immunodeficiency virus (HIV) or other infectious agents are absent, all products containing human source material should be considered potentially infectious and handled with the same precautions used with patient specimens.

This product contains <0.09% sodium azide as a preservative. Sodium azide may react with lead and copper plumbing to form potentially explosive compounds. Flush with excess water upon disposal.

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