Intended use

The Easicult S test is intended for detection and monitoring of anaerobic sulphide generating bacteria in various industrial environments.

Contents of the kit

Easicult S	Cat. No. 67687
Culture tubes	9 pcs
Capillary tubes	15 pcs
Labels	10 pcs
Instructions for use	1 pc

Typical formulation

Tryptone	Sodium lactate
Sodium sulphite	Sodium hydroxide
Iron II sulphate	Agar agar
Iron III citrate	Water
Magnesium sulphate	

Warnings and precautions

- Do not use the product beyond the expiry date marked on the kit.
- Handle the capillary only with forceps. Do not touch it with your hands or clothing and do not lay it down. Use a fresh capillary for each sample.
- Do not touch the unused growth media.
- Do not use the test if you notice discoloration of the growth medium or evidence of microbial growth.
- Because any growth on the Easicult S may be pathogenic, do not touch the growth.

Storage

Store Easicult S at room temperature (approx. 20°C/68°F) protected from draught, temperature fluctuations and light sources. Avoid storage near heat-generating appliances. Do not allow to freeze. The expiry date (year-month-date) is marked on the box.

Sampling and procedure (Fig. 1–6)

- Open the culture tube and stand it upright. Using forceps take one capillary from the container (if possible flame the jaws of the forceps before use). Caps should be placed on the table wide end down as illustrated.

- Dip the capillary into a tank, or
 - expose the capillary to a spray or running fluid, or
 - mix the sample in the container and dip the capillary filling it totally.
- Press the free end of capillary gently against a clean absorbent paper in order to empty the capillary.
- Insert the capillary into the culture medium (so that the capillary touches the bottom of the culture tube). Leave the capillary there.
- Close the culture tube firmly, fill in the label and affix it to the tube. Incubate the tube at about 35°C (95°F).
- Compare the degree of blackening in the tube with the model chart at regular, e.g. daily intervals for 5 days.

Interpretation of results

The growth of anaerobic sulfide generating bacteria in the Easicult S tube results in the formation of black iron sulphide. The blackening may begin at any point of the tube and depending on the degree of contamination eventually part of or the whole tube may become black. The result is not quantitative but shows the approximate level of contamination.

Once sulphide generation has commenced it will almost always increase fairly rapidly. Sulphide generation can be prevented by adequate aeration although aerobic micro-organisms will continue to flourish. Suppression of aerobic micro-organisms will usually suppress also sulfide generating organisms as the former are oxygen scavengers and create the oxygen free conditions needed by sulphide generating bacteria.

Disposal

Any growth on tests may be pathogenic. Used tests must therefore be disposed of according to local laws and regulations.

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Easicult® is a registered trademark of Orion Diagnostica Oy.