

EARLY DETECTION. RAPID DIAGNOSIS. BETTER OUTCOMES.

EARLY DETECTION. RAPID DIAGNOSIS. B

LY DETECTION. RAPID DIAGNOSIS. BETTER OUTCOMES.

EARLY DETECTION. RAPID DIAGNOSIS. BE

EARLY DETECTION. RAPID DIAGNOSIS. BETTER OUTCOMES.

EARLY DETECTION. RAPID DIAGNOSIS. BETTER OUTCOMES.

SIS. BETTER OUTCOMES. EARLY DETECTION. RAPID DIAGNOSIS

EARLY DETECTION. RAPID DIAGNOSIS. BETTER OUTCOMES.

DETECTION. RAPID DIAGNOSIS. BETTER OUTCOMES.

LY DETECTION. RAPID DIAGNOSIS. BETTER OUTCOMES.

DETECTION. RAPID DIAGNOSIS. BETTER O

EARLY DETECTION. RAPID DIAGNOSIS. BETTER OUTCOMES.

EARLY DETECTION. RAPID DIAGNOSIS

LY DETECTION. RAPID DIAGNOSIS. BETTER OUTCOMES.

TER OUTCOMES. EARLY DETECTION. RAPID DIAGNOSIS. BETTER

DIAGNOSIS. BETTER OUTCOMES.

EARLY DETECTION. RAPID DIAGNOSIS. BETTER OUTCOMES.

EARLY DETECTION. RAPID DIAGNOSIS. BETTER OUTCO

DIAGNOSIS. BETTER OUTCOMES.

CTION. RAPID DIAGNOSIS. BETTER OU

EARLY DETECTION. RAPID DIAGNOSIS. BETTER OUTCOMES.

TECTION. RAPID DIAGNOSIS. BETTER OUTCOMES.

EARLY DETECTION. R



EARLY DETECTION • RAPID DIAGNOSIS • BETTER OUTCOMES

Product Catalog

D DIAGNOSIS. BETTER OUTCOMES.

EARLY DETECTION. RAPID DIAGNOSIS. BETTER OUTCOMES.

EARLY DETECTION. RAPID DIAGNOSIS. BETTER

CTION. RAPID DIAGNOSIS. BETTER OUTCOMES.

CTION. RAPID DIAGNOSIS. BETTER OU

EARLY DETECTION. RAPID DIAGNOSIS. BETTER OUTCOMES.

EARLY DETECTION. RAPID DIAGNOSIS. BETTER OUTCOMES.

BETTER OUTCOMES. EARLY DETECTION. RAPID DIAGNOSIS. BET

EARLY DETECTION. RAPID DIAGNOSIS. BETTER OUTCOMES.

LifeSign provides high quality rapid point of care diagnostic products worldwide. We are proud to offer an extensive line of innovative products that are developed and manufactured in the USA under ISO, FDA, cGMP and CE mark guidelines.

Our product line includes Drugs of Abuse, Hygiene, Infectious Disease, Urology, Women's Health, Cardiac Care, GI Disorders and Industrial.

LifeSign's team of Sales, Technical and Customer Service professionals provide expert assistance and support.

CARDIAC CARE



LifeSign MI

For the Rapid Qualitative Detection of Cardiac Markers to Aid in the Diagnosis of Acute Myocardial Infarction in Human Whole Blood, Serum or Plasma.

- Qualitative CK-MB, Myoglobin, Troponin I
- Enables prompt Emergency Room decisions
- Easy to use procedure
- Results in 15 minutes
- Reduces unnecessary testing and admissions

ORDERING INFORMATION

Product	# OF TESTS	ITEM #
LifeSign MI Troponin I	20	60101
LifeSign MI Myoglobin/Troponin I	20	60301
LifeSign MI Myoglobin/CK-MB/Troponin I	20	60201
LifeSign MI Cardiac Control Set	2 X 1 ml	LSC-301

UROLOGY

Uricult®*

The Preferred In-Office Urine Culture System.

- Easy to perform; Dip - Incubate - Read Results
- Results within 18-24 hours (Incubation)
- Positive cultures can be sent to Reference Lab for additional studies if needed

*Trademark of Orion Diagnostica

ORDERING INFORMATION

Product	# OF TESTS	ITEM #
Uricult CLED/EMB	10	1000
Uricult CLED/MacConkey	10	1001
Uricult CLED/Polymixin/EMB	10	1002
Uricult CLED/Polymixin/MacConkey	10	1003



Status BTA

A Rapid Test for the Qualitative Detection of Bladder Tumor Associated Antigen in Human Urine.

- Indicated for use as an aid in monitoring bladder cancer patients in conjunction with cystoscopy
- **CLIA Waived**
- Accurate results with superior sensitivity & specificity
- Results in 5 minutes

ORDERING INFORMATION

Product	# OF TESTS	ITEM #
Status BTA	25	39025
Status BTA Control Set	2 X 1 ml	200019



UROLOGY



Uricult® Control Kit

Quality Control for Uricult Urine Culture Paddles. Each kit contains a positive and a negative control plus diluent.

ORDERING INFORMATION		
Product	# OF TESTS	ITEM #
Uricult Control Kit	2 x 99ml	200030

INFECTIOUS DISEASE

Status Flu A and B

For the Qualitative Detection of Influenza Type A and Type B Antigens Directly from Nasal Swab, Nasopharyngeal Swab and Nasal Aspirate/Wash Specimens.

- Flexible and flocked nasal swab for improved patient comfort and superior specimen collection
- Innovative flip test design with on board sample extraction
- Positive and Negative control swabs included
- Differentiates between Influenza A and B
- Individual pre-measured extraction reagent capsule for increased accuracy and ease of use



www.statusflu.com

ORDERING INFORMATION		
Product	# OF TESTS	ITEM #
Status Flu A and B	22	36022

Status Mono

For the Qualitative Detection of Infectious Mononucleosis Heterophile Antibodies in Whole Blood, Serum or Plasma.

- FDA cleared with no age restrictions
- **CLIA Waived** for Whole Blood
- Results within 8 minutes
- Requires only 25ul of whole blood or 10ul of serum/plasma
- Superior 99% sensitivity - specificity of 98.8%
- External positive and negative control included (item #68364 and 84M30)
- Safe and accurate plastic sample pipettes



ORDERING INFORMATION		
Product	# OF TESTS	ITEM #
Status Mono (WB/Serum/Plasma)	25	68364
Status Mono (WB/Serum/Plasma)	30	84M30
Status Mono (WB only)	30	84W30
Status Mono (WB only)	10	84W10
Status First Mono Control Set	2 X 0.5 ml	200014

INFECTIOUS DISEASE

Status Strep A Products

For the Rapid Detection of Group A Streptococcal Antigen
Directly from Throat Swab Specimens.

Status Strep A Flip

- **CLIA Waived**
- Sensitivity of 96.2%, specificity of 98.7%
- Results in 5-10 minutes
- External positive and negative controls included
- Innovative flip test device with on board sample extraction
- Individual, pre-measured sample extraction capsules



Status Strep A Strip

- **CLIA Waived**
- Sensitivity of 96.2%, specificity of 98.7%
- Results in 5-10 minutes
- External positive and negative controls included

Status Strep A Plus

- **CLIA Waived**
- Sensitivity of 96.2%, specificity of 98.7%
- Results in 5-10 minutes
- External positive and negative controls included
- Two convenient resealable canisters holding 25 test strips each



Status Accustrep Strip A

- Sensitivity of 96.2%, specificity of 98.7%
- Results in 5-10 minutes

ORDERING INFORMATION

Product	# OF TESTS	ITEM #
Status Strep A Flip (cassettes)	25	34125
Status Strep A Strip (pouch)	30	34130
Status Strep A Plus (cannister-strips)	50	34250
Status Accustrep Strip A (cassettes)	25	34025

INFECTIOUS DISEASE

Uni-Gold Recombigen HIV*

For Qualitative Detection of HIV-1/2 in Serum, Plasma and Whole Blood.

- **CLIA Waived** for whole blood
- Results in 10-12 minutes
- Sensitivity of 100%, Specificity of 99.7%

*Trademark of Trinity Biotech



ORDERING INFORMATION		
Product	# OF TESTS	ITEM #
Uni-Gold HIV	20	1206506
Uni-Gold HIV Control Set	2 X 0.5 ml	1206530

GI DISORDERS

Status H. pylori

For the Qualitative Detection of Anti-Helicobacter pylori Antibody in Human Whole Blood, Serum and Plasma.

- **CLIA Waived** for whole blood
- Requires only 25ul of whole blood or 10ul of serum or plasma
- Results in 10 minutes
- Sensitivity of 95.9% with overall agreement to biopsy histology



ORDERING INFORMATION		
Product	# OF TESTS	ITEM #
Status H. pylori	30	37030
Status H. pylori Control Set	2 X 1 ml	200013

Status iFOBT

Rapid Immunoassay for the Qualitative Detection of Immunochemical Fecal Occult Blood.
(Human Hemoglobin in Human Fecal Specimens)

- **CLIA Waived**
- No dietary or medication restrictions
- Increased patient compliance
- Results in 5 minutes
- Contains 30 test devices, 50 collection tubes, 50 mailing envelopes

www.statusifobt.com



ORDERING INFORMATION		
Product	# OF TESTS	ITEM #
Status iFOBT (complete kit)	30/50/50	38030-DTM
Status iFOBT Test Device Only	50	38050-D
Status iFOBT Sample Collection Tube Only	50	38050-T
Status iFOBT Mailing Envelope Only	50	38050-M
Status iFOBT Mailers & Collection Tubes	25	38025-TM
Status iFOBT Cassettes & Tubes	50	38050-DT
Status iFOBT Control Set	2 X 0.5 ml	200020

WOMEN'S HEALTH

Status hCG

For the Detection of hCG in Serum or Urine.

- **CLIA Waived** for urine
- Serum and urine 25mIU/ml sensitivity
- Results in 3-5 minutes
- Minimal susceptibility to the Hook effect by hCG β core fragment

Status hCG Strip

For the Detection of hCG in Urine.

- **CLIA Waived**
- Simple "dip and read" procedure
- 25mIU/ml sensitivity
- Results in 3-5 minutes
- Minimal susceptibility to the Hook effect by hCG β core fragment

StatusFirst hCG

For the Detection of hCG in Serum or Urine.

- Urine or Serum
- 25mIU/ml sensitivity
- Instrument read eliminates subjectivity (DxPress Reader)
- hCG Results in 5 minutes
- Minimal susceptibility to the Hook effect by hCG β core fragment



ORDERING INFORMATION

Product	# OF TESTS	ITEM #
Status hCG Urine/Serum	35	35035
Status hCG Urine	35	35135
Status hCG Strip Urine	35	35235
StatusFirst hCG (Read on DxPress Reader)	35	35335
Status hCG Control Set	2 X 5ml	135-100A

DxPRESS READER

DxPress® Reader

For the qualitative reading of LifeSign StatusFirst Test Devices - hCG and Drugs of Abuse.

- Eliminates visual read subjectivity for qualitative assays
- Easy to use
- Includes barcode scanner for test device and patient identification
- On board printer for hardcopy results
- RS-232 port for LIS connectivity
- Battery Power option
- Compact design



ORDERING INFORMATION

Product	# OF TESTS	ITEM #
DxPress® Reader (1 year warranty)	NA	LSR-2000
DxPress® Reader (3 year warranty)	NA	LSR-2003
StatusFirst hCG	35	35335
StatusFirst Drugs of Abuse (5 & 10 Panel)	25	30025
StatusFirst Troponin I (CE Marked - Export Only)	20	60501
StatusFirst D-Dimer (CE Marked - Export Only)	20	60601

DRUGS OF ABUSE

One Step Immunoassay for the Qualitative Detection of Drugs of Abuse and their Metabolites in Urine.

Status DS

- Simple one-step procedure
- Accurate results with 97-99% correlation to GC/MS
- Results within 5-10 minutes
- Flexible combination of drugs up to a 10 panel cassette



ORDERING INFORMATION

Product	# OF TESTS	ITEM #
Status DS TCA	10	20910
Status DS TCA	35	20935
Status DS 4 Panel THC/OPI/COC/MET	35	24135
Status DS 5 Panel THC/OPI/COC/MET/PCP	35	25135
Status DS 8 Panel THC/OPI/COC/AMP/BZO/ BAR/TCA/PCP	25	28025
Status DS 10 Panel MET/THC/OPI/COC/AMP/ BZO/BAR/TCA/PCP/MTD	10	21010
Status DS 10 Panel MET/THC/OPI/COC/AMP/ BZO/BAR/TCA/PCP/MTD	25	21025
Status DS Control Set	2 X 5 ml	21000



Status Cup Plus

- Split Specimen cup plus Status Stik Drug Test device
- Innovative cup designed for simple, accurate and secure test results
- 97-99% correlation to GC/MS
- Leak proof design eliminates contamination
- Three tamper evident seals ensure chain of custody

ORDERING INFORMATION

Product	# OF TESTS	ITEM #
Status Cup Plus/Status Stik THC/ OPI/COC/MET	35 cups/ 35 devices	74135
Status Cup Plus/Status Stik THC/ OPI/COC/MET/PCP	35 cups/ 35 devices	75135

Quick Screen 11 Drug Cup

- **CLIA Waived**
- Simple procedure with easy to read results
- Leak proof cup

ORDERING INFORMATION

Product	# OF TESTS	ITEM #
Quick Screen 11 MET/THC/OPI/COC/AMP/ BZO/BAR/TCA/PCP/MTD/MAD	25	60B02
Status DS Control Set	2 X 5 ml	21000



DRUGS OF ABUSE

Status Assure Cup & Drug Screen Card

- Simple, one-step procedure
- Results in 5 minutes
- Up to 11 drugs on one cup/card
- Adulterant strip included

ORDERING INFORMATION		
Product	# OF TESTS	ITEM #
Status Assure 6 Panel Cup THC/OPI-2000/ COC/MET/AMP/PCP	25	76025AA2K (cup)
Status Assure 11 Panel Cup THC/OPI-300/ OXY/MDMA/TCA/BZO/BAR/COC/MET/AMP/ PCP	25	71125AA (cup)
Status Assure 11 Panel Cup THC/OPI-2000/ OXY/MDMA/TCA/BZO/BAR/COC/MET/AMP/ PCP	25	71125AA2K (cup)
Status Assure 11 Panel Card THC/OPI-300/ OXY/MDMA/TCA/BZO/BAR/COC/MET/AMP/ PCP	25	11125AA (card)
Status Assure 11 Panel Card THC/OPI-2000/ OXY/MDMA/TCA/BZO/BAR/COC/MET/AMP/ PCP	25	11125AA2K (card)
Status DS Control Set	2 X 5 ml	21000



StatusFirst DOA

- Instrument read eliminates subjectivity
- Accurate qualitative results in 5 minutes
- Available in 5 and 10 drug panels
- Results stored and printed for documentation

Instrument Read-Qualitative DXpress



ORDERING INFORMATION		
Product	# OF TESTS	ITEM #
Status First DOA 10 MET/THC/OPI/COC/AMP/ BZO/BAR/TCA/PCP/MTD*	25	30025
Status First DOA 5 MET/THC/OPI/COC/PCP*	25	30525
Status DS Control	2 X 5 ml	21000

**Read on DXpress® Reader*

ALCOHOL

Q.E.D. A150

For the Rapid Quantitative Determination of Alcohol in Saliva.

- **CLIA Waived**
- Simple two step procedure
- Accurate results; built-in QA spot™ assures quality control

ORDERING INFORMATION		
Product	# OF TESTS	ITEM #
Q.E.D.A150	10	A-150



NICOTINE

Status Nicotine

Rapid Immunoassay for the Detection of Cotinine, Nicotine Metabolite in Urine.

- Simple one step procedure
- Results in 10 minutes
- 500ng/ml cut-off



ORDERING INFORMATION		
Product	# OF TESTS	ITEM #
Status DS Nicotine	35	21735
Status DS Controls	2 x 5 ml	21000

HYGIENE

Hygicult™*

Rapid Hygiene Culture System.

- Ready to use for on site hygiene monitoring
- Store at room temperature
- Hinged slide for convenient surface sampling
- AOAC approved
- Incubate at 36°C for 24 hours



ORDERING INFORMATION		
Product	# OF TESTS	ITEM #
Hygicult™ TPC	10	68010
Hygicult™ Y&F	10	68013
Hygicult™ E/Beta-Gur	10	68267

*Trademark of Orion Diagnostica

Clean Card® Pro*

Detects Protein Residues on Surfaces.

- Easy to use surface hygiene test pad
- Results in 30 seconds
- Long shelf life; stored at room temperature
- No instrument or incubation required

ORDERING INFORMATION		
Product	# OF TESTS	ITEM #
Clean Card® Pro 25 (includes water spray bottle)	25	133148
Clean Card® Pro 50	50	133149



INDUSTRIAL

Easicult™*

Convenient Industrial Culture System Identifying Bacteria & Fungi in Oils, Water, Fuels and Other Applications.

- Ready to use on site test for monitoring microbial contamination
- Easy result interpretation
- Decreases costly disposal of unusable fluids
- Incubate at 36°C for up to 24 hours
- Option to incubate at room temperature for up to 5 days



ORDERING INFORMATION

Product	# OF TESTS	ITEM #
Easicult® TTC	10	4030
Easicult® Combi	10	4060

*Trademark of Orion Diagnostica

INCUBATOR

Cultura® M Incubator

Small, Durable and Reliable Incubator. Suitable for the Incubation of Dip Paddles.

- Lightweight and portable
- Temperature range from 25°C to 48°C
- Holds 18 dip paddles



ORDERING INFORMATION

Product	# OF TESTS	ITEM #
Cultura® M Incubator	NA	05866

lifeSign

EARLY DETECTION • RAPID DIAGNOSIS • BETTER OUTCOMES



lifeSign

A PBM Group Company

85 Orchard Road • Skillman • New Jersey • 08558

800.526.2125 • 732.246.3366 • Fax: 732.246.0570

www.lifesignmed.com