




lifeSign
THE DXPRESS READER





DXpress® Reader

A multifunctional reader with user friendly features to deliver fast and accurate results in the laboratory or at the point-of-care

- Easy to use; includes a barcode scanner for easy patient I.D. input
 - LOT specific test device information on data chip
 - Operator & QC lockout features
 - QC controls built into the reader software & test device
 - Stores & manages up to 1,000 patient and 60 QC results
-
- On-board printer for hardcopy results
 - RS-232 port for LIS connectivity when downloading patient results
 - Instrument read eliminates visual subjectivity
 - Lightweight for portability, less than 3 lbs
 - Single or batch testing with on-board timer



ORDERING INFORMATION

PRODUCT	ITEM #
DXpress® Reader	LSR-2000 (1 year warranty)
DXpress® Reader	LSR-2003 (3 year warranty)

StatusFirst® CHF NT-proBNP

For the rapid quantitative detection of NT-proBNP aiding in the diagnosis of congestive heart failure



- Quantitative results in 15 minutes with DXpress® Reader
- Easy to use procedure; Add 3 drops of specimen into device
- EDTA plasma sample
- Measures a range of 20 to 5000 pg/mL
- Analytical sensitivity of 20 pg/mL


ORDERING INFORMATION

PRODUCT	# OF TESTS	ITEM #
StatusFirst® CHF NT-proBNP	20	60401
CONTROL PRODUCT	# OF TESTS	ITEM #
StatusFirst® CHF NT-proBNP Controls	2 x 1 mL	LSCHF-300C

StatusFirst® DOA

One step immunoassay for the qualitative detection of drugs of abuse and their metabolites in urine

- FDA cleared
- Simple one-step procedure
- Instrument read eliminates subjectivity
- Accurate results in 5 minutes
- Available in 5 & 10 drug panels
- Results stored and printed for documentation
- Walk away procedure accommodates busy labs and workflow interruptions



ORDERING INFORMATION

PRODUCT	# OF TESTS	ITEM #
StatusFirst® DOA 5 Panel THC/COC/OPI/MET/PCP	25	30525
StatusFirst® DOA 10 Panel MET/OPI/COC/THC/PCP/BZO/BAR/MTD/TCA/AMP	25	30025
CONTROL PRODUCT	# OF TESTS	ITEM #
StatusFirst® DOA 10 Controls	2 x 5 mL	21000

StatusFirst® hCG

For the detection of hCG in serum or urine



- Urine or Serum
- Superior sensitivity
- Simple procedure
- Instrument read eliminates subjectivity
- hCG results in 3-5 minutes

ORDERING INFORMATION

PRODUCT	# OF TESTS	ITEM #
StatusFirst® hCG Urine/Serum	35	35335
CONTROL PRODUCT	# OF TESTS	ITEM #
StatusFirst® hCG Urine Controls	2 x 5 mL	135-100A

StatusFirst® TnI Troponin I*

For the rapid quantitative determination of Troponin I in human whole blood, plasma or serum

- Whole Blood, Plasma or Serum
- Easy to use - 130 µL of specimen into test device
- Quantitative results in 15 minutes in DXpress® Reader



ORDERING INFORMATION

PRODUCT	# OF TESTS	ITEM #
StatusFirst® TnI	20	60501
CONTROL PRODUCT	# OF TESTS	ITEM #
StatusFirst® MI Controls	2 x 3 mL	LSC-300C

StatusFirst® D-Dimer*

For the quantitative determination of D-Dimer as an aid in the diagnosis of Deep Vein Thrombosis (DVT), Pulmonary Embolism (PE) and Disseminated Intravascular Coagulation (DIC)



- Whole Blood or Plasma
- Easy to use - 130 µL of specimen into text device
- Quantitative results in 15 minutes in DXpress® Reader
- Measures a range of 60 to 5000 ng/mL (.06 mg/L - 5 mg/L)
- Sensitivity 100% at cut off 250 ng/mL (0.25 mg/L)
- Negative predictive value 100%

ORDERING INFORMATION

PRODUCT	# OF TESTS	ITEM #
StatusFirst® D-Dimer	20	60601
CONTROL PRODUCT	# OF TESTS	ITEM #
StatusFirst® CHF NT-proBNP Controls	2 x 1 mL	LSCHF-300C





lifeSign

EARLY DETECTION • RAPID DIAGNOSIS • BETTER OUTCOMES



**For more information or to
place an order, contact:**



lifeSign

A PBM Group Company

71 Veronica Avenue • Somerset • New Jersey • 08873

800.526.2125 • 732.246.3366 • Fax: 732.246.0570

www.lifesignmed.com