

Status

# hCG

lifeSign

Rapid immunoassay for the qualitative detection of human chorionic gonadotropin

## Status hCG Cassette

For the detection of hCG in serum or urine

*The one to choose for exceptional sensitivity and early pregnancy detection*

### Features

Superior sensitivity, 25mIU/mL, regardless of specimen for early pregnancy detection

Clear, definitive result lines are easy to read

Simple one-step procedure, just add 3 drops of urine or serum and read in 5 minutes

CLIA Waived for urine



Status hCG

hCG Strip

## Status hCG Strip

For the detection of hCG in urine

*The most economical choice*

### Features

"Dip and read" simplicity makes it easy to use

25 mIU/mL sensitivity for early pregnancy detection

Compact strip size, reduces biohazard disposal

Results in minutes

CLIA Waived

Other **Status** Testing  
Solutions available from LifeSign:

Strep A

Cardiac Markers

H. pylori

Drugs of Abuse

Mono

For more information about  
Status hCG,  
or to place an order, contact:

### ORDERING INFORMATION

PRODUCT	SENSITIVITY	KIT SIZE	ITEM #	CPT #
Status hCG urine/serum	25 mIU/mL serum 25 mIU/mL urine	35 tests	35035	84703 81025
Status hCG urine only	25 mIU/mL urine	35 tests	35135	81025QW
Status hCG Strip urine only	25 mIU/mL urine	35 tests	35235	81025QW

lifeSign

A PBM Group Company

71 Veronica Avenue  
Somerset, New Jersey 08873  
800.526.2125 732.246.3366  
Fax: 732.246.0570  
www.lifesignedmed.com

LifeSign is a registered trademark of Princeton BioMeditech Corporation. Printed in USA. All rights reserved.