

DRUGS OF ABUSE

Rapid immunoassay for the qualitative detection of Methamphetamine, Opiates, Cocaine, THC (Cannibus), Phencyclidine, Benzodiazepines, Barbiturates, Methadone, Tricyclic Antidepressants, Amphetamine, and/or their Metabolites in Urine

StatusFirst® DOA

For the detection of DOA in urine

Features:

- Simple procedure, 2 drops of urine
- Instrument read eliminates subjectivity
- Accurate results in 5 minutes
- FDA Cleared



StatusFirst® DOA

DXpress® Reader

A multifunctional reader which delivers accurate and efficient results in point-of-care situations

Features:

- Easy to use; includes a barcode scanner
- Stores & manages up to 1,000 patient results
- On-board printer for hardcopy results
- RS-232 port for LIS connectivity for downloading patient results
- Lightweight for portability, less than 3 lbs



DXpress® Reader

ORDERING INFORMATION

PRODUCT	KIT SIZE	ITEM #	CPT CODE
StatusFirst® DOA 5 Panel (MET/OPI/COC/THC/PCP)	25 tests	30525	80101
StatusFirst® DOA 10 Panel (MET/OPI/COC/THC/PCP/BZO/BAR/MTD/TCA/AMP)	25 tests	30025	80101
DXpress® Reader (1 year warranty)		LSR-2000	

For more information about StatusFirst® DOA, DXpress® Reader or to place an order, contact:

lifeSign

A PBM Group Company
 71 Veronica Avenue
 Somerset, New Jersey 08873
 800.526.2125 732.246.3366
 Fax: 732.246.0570
www.lifesignmed.com

Currently Available Products:

- NT-proBNP - Quantitative
- hCG - Qualitative
- Drugs of Abuse - Qualitative

Products Coming Soon:

- Troponin I
- D-Dimer
- Infectious Disease

DXpress is a registered trademark of LifeSign. LifeSign is a registered trademark of Princeton BioMeditech Corporation. Printed in USA. All rights reserved.