

EARLY DETECTION. RAPID DIAGNOSIS. BETTER OUTCOMES.

EARLY DETECTION. RAPID DIAGNOSIS. BE

EARLY DETECTION. RAPID DIAGNOSIS. BETTER OUTCOMES.

EARLY DETECTION. RAPID DIAGNOSIS. BET



# lifeSign

## PRODUCT CATALOG

EARLY DETECTION. RAPID DIAGNOSIS. BETTER OUTCOMES.

EARLY DETECTION. RAPID DIAGNOSIS. BETTER OUTCOMES.

EARLY DETECTION. RAPID DIAGNOSIS.

BETTER OUTCOMES. EARLY DETECTION. RAPID DIAGNOSIS. BETTER

DIAGNOSIS. BETTER OUTCOMES.

EARLY DETECTION. RAPID DIAGNOSIS. BETTER OUTCOMES.

EARLY DETECTION. RAPID DIAGNOSIS. BETTER OUTCO

DIAGNOSIS. BETTER OUTCOMES.

CTION. RAPID DIAGNOSIS. BETTER OUT

EARLY DETECTION. RAPID DIAGNOSIS. BETTER OUTCOMES.

TECTION. RAPID DIAGNOSIS. BETTER OUTCOMES.

EARLY DETECTION. RA

BETTER OUTCOMES. EARLY DETECTION. RAPID DIAGNOSIS. BETTER

CTION. RAPID DIAGNOSIS. BETTER OUTCOMES.

EARLY DETECTION. RAPID DIAGNOSIS. BETTER OUTCO

DIAGNOSIS. BETTER OUTCOMES.

CTION. RAPID DIAGNOSIS. BETTER OUT

**RAPID DIAGNOSIS • EARLY DETECTION • BETTER OUTCOMES**

CTION. RAPID DIAGNOSIS. BETTER OUT

EARLY DETECTION. RAPID DIAGNOSIS. BETTER OUTCOMES.

EARLY DETECTION. RAPID DIAGNOSIS. BETTER OUTCOMES.

BETTER OUTCOMES. EARLY DETECTION. RAPID DIAGNOSIS. BETT

EARLY DETECTION. RAPID DIAGNOSIS. BETTER OUTCOMES.

EARLY DETECTION. RAPID DIAGNOSIS. BETT



## Status™ COVID-19/Flu A&B

*A Rapid Immunoassay for the Simultaneous Direct Detection and Differential Diagnosis of SAR-CoV-2, Influenza Type A and Type B Antigen from anterior nasal and nasopharyngeal swab specimens*



- COVID-19 Anterior nasal swab specimen – Sensitivity 93.8%, Specificity 100%
- COVID-19 Nasopharyngeal – Sensitivity 93.1%, Specificity 100%
- Flu A – Sensitivity 91.4%, Specificity 95.7%
- Flu B – Sensitivity 87.6%, Specificity 95.9%
- FDA Emergency Use Authorization (EUA)
- Visually read in 15 minutes
- Flocked nasopharyngeal swab for superior specimen collection and patient comfort

ORDERING INFORMATION		
PRODUCT	#OF TESTS	ITEM#
Status™ Covid-19/Flu A&B	25	33225

## Status™ Flu A and B

For the Qualitative Detection of Influenza Type A and Type B Antigens Directly from Nasal Swab, Nasopharyngeal Swab and Nasal Aspirate/Wash Specimens.

- **CLIA Wavied** for Swab Specimens
- Flexible and flocked nasal swab for improved patient comfort and superior specimen collection
- Innovative flip test design with no board sample extraction
- Positive and Negative control swabs included
- Differentiates between Influenza A and B
- Individual pre-measured exaction reagent capsule for increased accuracy and ease of use



Flocked Swab

### ORDERING INFORMATION

PRODUCT	#OF TESTS	ITEM#
Status™ Flu A&B	25	36025
Nasal Aspirate kit (44 pipettes, 1 reagent)	22	510AS
Status Flu A&B Controls	5 sets	36022CS
Status A&B nasal swab refill	52	36052

## Status™ Strep A Products

For the detection of Group A Streptococcal Antigen Directly from Throat Swabs Specimens

### Status™ Strep A Flip

- **CLIA Wavied**
- Sensitivity of 96.2%, specificity of 98.7%
- Results in 5-10 minutes
- External positive and negative controls included
- Innovative flip test device with on board sample extraction
- Individual, pre-measured sample extraction capsules



### Status™ Strep A Strip

- **CLIA Wavied**
- Sensitivity of 96.2%, specificity of 98.7%
- Results in 5-10 minutes
- External positive and negative controls included



## Status™ Strep A Plus

- **CLIA Wavied**
- Sensitivity of 96.2%, specificity of 98.7%
- Results in 5-10 minutes
- External positive and negative controls included
- Two convenient resealable canisters holding 25 test strips each



ORDERING INFORMATION		
PRODUCT	#OF TESTS	ITEM#
Status Strep A Flip (cassettes)	25	34125
Status Strep A Strip (pouch)	30	34130
Status Strep A Plus (canister-strips)	50	34250

## Status™ Mono

For the Qualitative Detection of Infectious Mononucleosis Heterophile Antibodies in Whole Blood, Serum or Plasma.

- FDA cleared with no age restrictions
- **CLIA Wavied** for Whole Blood
- Results in 8 minutes
- Requires only 25ul of whole blood or 10 ul serum/ plasma
- Superior 99%, sensitivity - specificity of 98.8%
- External positive and negative control included (item #68364 and 84M30)
- Safe and accurate plastic sample pipettes



ORDERING INFORMATION		
PRODUCT	#OF TESTS	ITEM#
Status Mono (WB/Serum/Plasma)	25	68364
Status Mono (WB/Serum/Plasma)	30	84M30
Status Mono (WB only)	30	84W30
Status Mono (WB only)	10	84W10
Status Mono Control Set	2 x 1 ml	200014

## Status™ H. pylori

For the Qualitative Detection of Anti-Helicobacter pylori Antibody in Human Whole Blood, Serum and Plasma.

- **CLIA Wavied** whole blood
- Requires only 25ul of whole blood or 10 ul of serum or plasma
- Results in 10 minutes
- Sensitivity of 95.9% with overall agreement to biopsy histology



### ORDERING INFORMATION

PRODUCT	#OF TESTS	ITEM#
Status™ H.pylori	30	37030
Status™ H.pylori Control Set	2 X 1 ml	200013

## Status™ iFOBT

Rapid Immunoassay for the Qualitative Detection of Immunochemical Fecal Occult Blood. (Human Hemoglobin in Human Fecal Specimens)

- **CLIA Wavied**
- No dietary or medication restrictions
- Increased patient compliance
- Results in 5 minutes
- Contains 30 test devices, 50 collection tubes, 50 mailing envelopes



### ORDERING INFORMATION

PRODUCT	#OF TESTS	ITEM#
Status iFOBT (complete kit)	30/50/50	38030 -DTM
Status iFOBT Sample Collection Tube Only	50	38050 -T
Status iFOBT Mailing Envelope Only	50	38050-M
Status iFOBT Mailer & Collection Tubes	25	38025-TM
Status iFOBT Device & Tubes	50	38050-DT
Status iFOBT Control Set	2 X 1 ml	200020

## Uricult®

*The Preferred In-Office Urine Culture System.*

- Easy to perform; Dip - Incubate - Read Results
- Results within 18-24 hours (Incubation)
- Positive culture can be sent to Reference Lab for additional studies if needed

• Trademark of Orion Diagnostica



### ORDERING INFORMATION

PRODUCT	#OF TESTS	ITEM#
Uricult® CLED/EMB	10	1000
Uricult® CLED/Polymixin/EMB	10	1002
Uricult® CLED/Polymixin/MacConkey	10	1003



## Uricult® Control Kit

Quality Control for Uricult Urine Culture Paddles. Each kit contains a positive and a negative control plus diluent.

### ORDERING INFORMATION

PRODUCT	#OF TESTS	ITEM#
Uricult® Control Kit	2 x 99ml	200030

## Status™ BTA

*A Rapid Test for the Qualitative Detection of Bladder Tumor Associated Antigen in Human Urine.*

- **CLIA Wavied**
- Indicated for use as an aid in monitoring bladder cancer patients in conjunction with cystoscopy
- Accurate results with superior sensitivity & specificity
- Results in 5 minutes



### ORDERING INFORMATION

PRODUCT	#OF TESTS	ITEM#
Status BTA	25	39025
Status BTA Control Kit	2 X 2 ml	200019

## CARDIAC CARE



### LifeSign™ MI

For the Rapid Qualitative Detection of Cardiac Markers to Aid in the Diagnosis of Acute Myocardial Infarction in Human Whole Blood, Serum or Plasma.

- Qualitative CK-MB, Myoglobin, Troponin
- Enables prompt Emergency Room decisions
- Easy to use procedure
- Result in 15 minutes
- Reduces unnecessary testing and admissions

#### ORDERING INFORMATION

PRODUCT	#OF TESTS	ITEM#
LifeSign MI Troponin I	20	60101
LifeSign MI Myoglobin (CK-MB)/Troponin I	20	60201
LifeSign MI Cardiac Control Kit	2 X 1 ml	LSC-301(a)

## WOMEN'S HEALTH



### Status™ hCG

For the Detection of hCG in Serum or Urine.

- **CLIA Waved** for Urine
- Serum and urine 25mIU/ml sensitivity
- Results in 3-5 minutes
- Minimal susceptibility to the Hook effect hCG  $\beta$  core fragment

### Status™ hCG Strip

For the Detection of hCG in Urine.

- **CLIA Waved**
- Simple "dip and read" procedure
- 25mIU/ml sensitivity
- Results in 3-5 minutes
- Minimal susceptibility to the Hook effect hCG  $\beta$  core fragment



#### ORDERING INFORMATION

PRODUCT	#OF TESTS	ITEM#
Status hCG Urine/Serum (cassette)	35	35035
Status Urine (cassette)	35	35135
Status hCG Strip	35	35235
Status hCG Control Set	2 X 2.5 ml	135-100A




## DRUGS OF ABUSE

### Status™ DS

One Step Immunoassay for the Qualitative Detection of Drugs of Abuse and their Metabolites in Urine.

- Simple one-step procedure
- Accurate results with 97-99% correlation to GC/MS
- Results in 5-10 minutes



ORDERING INFORMATION		
PRODUCT	#OF TESTS	ITEM#
Status DS TCA	35	20935
Status DS Control Set	2 X 5 ml	21000
Status DS 10 panel MET/THC/OPI/COC/AMP/BZO/ BAR/TCA/PCP/MTD	25	20125

## NICOTINE

### Status™ Nicotine

Rapid Immunoassay for the Detection of Cotinine, Nicotine Metabolites in Urine.

- Simple one-step procedure
- Accurate results with 97-99% correlation to GC/MS
- Results in 10 minutes
- 500mg/ml cut-off



ORDERING INFORMATION		
PRODUCT	#OF TESTS	ITEM#
Status DS Nicotine	35	21735
Status DS Control Set	2 X 5 ml	21000

**Hygicult®**

*Rapid Hygiene Culture System.*

- Ready to use for on site hygiene monitoring
- Store at room temperature
- Hinged slide for convenient surface sampling
- AOAC approved
- Incubate at 36°C for 24 hours



**ORDERING INFORMATION**

PRODUCT	#OF TESTS	ITEM#
Hygicult™ TPC	10	68010
Hygicult™ Y&F	10	68013

**Easicult®**

*Convenient Industrial Culture System Identifying Bacteria & Fungi in Oils, Water, Fuels and Other Applications.*

- Ready to use on site test for Monitoring microbial contamination
- Easy result interpretation
- Decease costly disposal of unsable fluids
- Incubate at 36oC for up to 24 hours
- Option to incubate at room tempreature for up to 5 days



**ORDERING INFORMATION**

PRODUCT	#OF TESTS	ITEM#
Easicult™ TTC	10	4030
Easicult™ Combi	10	4060

# INCUBATOR



## Cultura® M Incubator

*Small, Durable and Reliable Incubator. Suitable for the Incubation of Dip Paddles.*

- Lightweight and portable
- Temperature range from 25°C to 48°C
- Holds 18 dip paddles

### ORDERING INFORMATION

PRODUCT	#OF TESTS	ITEM#
Cultura®M Incubator	NA	05866

# lifeSign

EARLY DETECTION • RAPID DIAGNOSIS • BETTER OUTCOMES



lifeSign

**A PBM Group Company**

85 Orchard Road • Skillman • New Jersey • 08558

800.526.2125 • 732.246.3366 • Fax: 732.246.0570

[www.lifesignmed.com](http://www.lifesignmed.com)